

Freestyle Libre (Interstitial Glucose) Monitor:

Your pet has a Freestyle Libre interstitial glucose monitoring system in place. This system monitors interstitial glucose continually. To get an initial reading, scan the reader (or your smartphone using the free app) over the sensor placed on your pet's shoulder. Please take one reading every 8 hours in order to avoid gaps in data collection (the device can only "reach back" 8 hours at a time). Each sensor lasts up to 14 days.

Please set up an online Libreview account (libreview.com) and share your account with saim@tier1vet.com. Send an email update to saim@tier1vet.com including appetite, water intake, urination habits, and behavior after the sensor has been in place for 7-10 days. This email will allow us to interpret your pets diabetic control in light of clinical signs at home.

Please set up a patient account at:

<https://www.libreview.com/>

Free app download (Freestyle libre 2 app) for at home monitoring:

<https://www.freestyle.abbott/us-en/products/freestyle-libre-app.html>

FreeStyle Libre 2 app



FreeStyle Libre 2 app*

